




Automatic Statistical Process Control Analysis of Pills, Capsules & Tablets



Package Order# LA230P-PF

LA230P : 60/120/230g x 0.1/0.2/0.5mg
YTF02-115NA PILL FEEDER

Includes:

- All necessary cables
- SPC software
- Stainless steel draft shield door (left hand side)
- Weighing bowl
- YDP20-OCE statistical printer with T-Connector

List Price- \$12,855.00 USD

Features & Benefits:

- Simple SPC program produces sophisticated statistics according to USP guidelines
- Pill feeder offers variable vibration settings for any type of pill or capsule
- Plug & play setup allows user to begin analysis immediately
- Eliminates the need for manual weighing of pills and the risk of inherent errors in data recording
- The supplied printer allows configurable printout of results (can be connected to PC)

Sartorius Mechatronics Corporation
5 Orville Drive
Bohemia, NY 11716

Phone: 800.635.2906
Fax: 631.254.4252

Sartorius Mechatronics Canada Inc.
2179 Dunwin Dr. Units #4/5
Mississauga, Ontario L5L 1X2

Phone: 800.668.4234
Fax: 905.569.7021

www.sartorius-mechatronics.com